1. Indication
audioprint® HPR-1 from pro3dure is a methacrylate based resin for 3D-printing systems (SLA, DLP) for the production of customized hearing protection.

2. Contraindication
audioprint® HPR-1 is contraindicated for every application that is not part of the indication (see above).

3. Requirements
Software – Information available from:
Cyfex AG · Binzmühlestrasse 15, 8050 Zürich, Switzerland
3SHAPE A/S · Holmens Kanal 7, 1060 Copenhagen, Denmark

Hardware (3D Printing) – Information available from:
ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

Hardware (Post Curing) – Information available from:
pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany

Additional specifications on www.pro3dure.com.

4. Material
Contains: Alkoxylated bisphenole-A-dimethacrylate, initiators, dyes, stabilisers, and pigments.

5. Geometric presetting
According to the appropriate design specifications (templates).

6. Material parameters
Radiation penetration depth controllable by exposure time.

7. Manufacturing process (fig. 1–10, page 1)

1. Prepare data (CAD & build preparation).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (with IPA ≥ 97 % or equivalent cleaning agent) approx. 5 min. in an ultrasonic bath or equivalent device – precleaning recommended.
8. Dry parts, until there are no residues of IPA or equivalent cleaning agent.
9. Post curing (10 min.): inert atmosphere recommended (use adequate light curing device).
10. Finish parts.

8. Finishing processes
Polishing / Lacquering audioprint® L2-FLX

9. Notice
Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.

According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

Caution: Polymerized resins are chemically resistant – avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

10. Hazard statements
Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Harmful to aquatic life with long lasting effects.

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